mikrobiologie labor technik

Hektoen Agar

PRINCIPIU

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
Hektoen Agar

PRINCIPLE

Hektoen Agar is used to isolate and differentiate pathogen Enterobacteria. The inhibitor effect of bile salts associated to excellent nutrient properties makes it an ideal medium to detect Salmonella and Shigella. Proteus do not spread on this medium. The principle of this medium relies on the fermentation of three carbohydrates : salicine, lactose and saccharose present in the medium.

FORMULA

Components	g/l
Meat peptone	12.00
Yeast extract	3.00
Bile salts	9.00
Lactose	12.00
Saccharose	12.00
Salicine	2.00
Sodium chloride	5.00
Sodium thiosulfate	5.00
Ferric citrate	1.50
Bromothymol blue	0.064
Acid Fuchsine	0.10
Agar	15.00
Final pH : 7.6 ± 0.2 at 25°C	

PREPARATION

Suspend 76.6 grams of powder in one litre of purified water. Bring to the boil until complete dissolution, and boil for a few minutes. DO NOT AUTOCLAVE.

PROCEDURE

Liquefy the medium in a boiling water bath then cool to 44-47°C. Pour into sterile

Petri plates, let them set on a flat surface. Dry in an incubator with the lids slightly open. Inoculate the plate by streaks from enrichment media used for the detection of Enterobacteria. Incubate for 24 and 48 hours at 37°C.

RESULTS

Microorganisms that ferment one of the 3 carbohydrates will grow as salmon pink colonies, the others grow as blue or green colonies. In presence of sodium thiosulfate, H₂ S producing microorganisms will give in presence of ferric citrate a black centre to their colonies.

Colonies characteristics :

Salmon / yellow	<i>Arizona, Citrobacter, Enterobacter, Escherichia, Klebsiella, Serratia, Y. enterocolitica.</i>
Black centred Salmon / yellow	<i>Proteus vulgaris.</i>
Black centred green colonies	<i>Proteus mirabilis, Salmonella.</i>
Green or bluish colonies	<i>Morganella morganii, Proteus rettgeri, Providencia, Salmonella H2S -, Shigella.</i>
Small blue or brownish colonies	<i>Pseudomonas.</i>

LIMITS AND PRECAUTIONS

Some strains of Vibrio can grow on this medium as salmon yellow colonies.

BIBLIOGRAPHY

1. King S. and Metzger W.I. 1968. A new plating medium for the isolation of enteric pathogens, II. Comparison of Hektoën Enteric Agar with SS and EMB Agar. Appl. Microbiol, 16:579-581.
2. Taylor W.I. and Schelhaut D. 1971. Comparison of Xylose Lysine Desoxycholate Agar, Hektoën Enteric Agar, Salmonella-Shigella Agar, and Eosin Methylene Blue Agar with stool specimens. Appl. Microbiol. 21:32- 37.
3. King S. and Metzger W.I. 1968. A new plating medium for the isolation of

enteric pathogens. I. Hektoën Enteric Agar. Appl. Microbiol. 16:577-578.

PACKAGING

Dehydrated

(To store between 1 and 30°C)

140132A: Flask of 500 g

140132L: Flask of 500 g

Ready to use medium

(To store between 2 and 25°C)

130132: Pack of 3 flasks of 200 ml

110132: Pack of 20 plates 90 mm Ø